

Anti-Cab39 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y376

Item Number: CM0018813



Introduction

High-quality anti-CAB39 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Targets calcium-binding protein 39 (MO25alpha), a key stabilizer of the STK11/LKB1 complex.

[Learn More](#)

Parameter	Value
Product Name	CAB39 Rabbit pAb
Remarks/Alias	MO25; CGI-66; CAB39
Species	Human
Gene ID (Human)	51719
Gene ID	51719
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-341 of human CAB39 (NP_057373.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y376
Protein Weight	40kDa
Shipping	Ice bag