


Anti-Cab39L Rabbit Polyclonal Antibody For Wb, Elisa - Q9H9S4

Item Number: CM0029148



Introduction

CAB39L rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse CAB39L (39kDa). Targets calcium-binding protein 39-like (MO25beta), a STK11/LKB1 activator complex component.

[Learn More](#)

Parameter	Value
Product Name	CAB39L Rabbit pAb
Remarks/Alias	MO2L; MO25-BETA; bA103J18.3; CAB39L
Species	Human
GeneID (Human)	81617
GeneID	81617
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-337 of human CAB39L (NP_112187.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H9S4
Protein Weight	39kDa
Shipping	Ice bag