

Anti-Cables1 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tdn4

Item Number: CM0016774

Introduction

CABLES1 Rabbit pAb for WB, ELISA. Validated in human, mouse. CDK5/ABL1 substrate 1, key in neuronal development and cell cycle. From CamelBio, one-stop IVD raw materials supplier.

[Learn More](#)

Parameter	Value
Product Name	CABLES1 Rabbit pAb
Remarks/Alias	CABL1; IK3-1; CABLES; HsT2563; CABLES1
Species	Human
GeneID (human)	91768
GeneID	91768
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-368 of human CABLES1 (NP_612384.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TDN4
Protein Weight	68kDa
Shipping	Ice bag