


Anti-Cabp1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nzu7

Item Number: CM0020611



Introduction

High-quality rabbit polyclonal antibody against human CABP1 (Calbrain/Caldendrin), validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat CABP1 with a predicted molecular weight of 40kDa.

[Learn More](#)

Parameter	Details
Product Name	CABP1 Rabbit pAb
Remarks/Alias	CALBRAIN; HCALB_BR; CABP1
Species	Human
GeneID	9478
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CABP1 (NP_001028849.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZU7
Protein Weight	40kDa
Shipping	Ice bag