


## Anti-Cabp7 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86V35

Item Number: CM0013851



### Introduction

Rabbit polyclonal antibody against CABP7 (Calneuron II) validated for Western blot, IHC-P, and ELISA in human, mouse, and rat. Targets a Golgi trafficking regulator. Ideal for cell biology studies.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | CABP7 Rabbit pAb                                                                                                    |
| Remarks/Alias    | CALN2; CABP7                                                                                                        |
| Species          | Human                                                                                                               |
| GeneID (human)   | 164633                                                                                                              |
| GeneID           | 164633                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 100-180 of human CABP7 (NP_872333.1). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q86V35                                                                                                              |
| Protein Weight   | 24kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |