


Anti-Cacna1D Rabbit Polyclonal Antibody For Wb And Elisa - Q01668

Item Number: CM0019544



Introduction

Rabbit polyclonal antibody against human CACNA1D (Cav1.3), reactive with human and rat, suitable for Western blot and ELISA. Explore L-type calcium channel function in excitable cells.

[Learn More](#)

Parameter	Value
Product Name	CACNA1D Rabbit pAb
Remarks/Alias	CACH3; CACN4; PASNA; SANDD; Cav1.3; CCHL1A2; CACNL1A2; CACNA1D
Species	Human
Gene ID (Human)	776
Gene ID	776
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1882-2181 of human CACNA1D (NP_000711.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q01668
Protein Weight	245kDa
Shipping	Ice bag