

Anti-Cacna1E Rabbit Polyclonal Antibody For Wb, Elisa - Q15878

Item Number: CM0015021



Introduction

CACNA1E Rabbit Polyclonal Antibody validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Raised against recombinant protein, suitable for R-type calcium channel research.

[Learn More](#)

Parameter	Value
Product Name	CACNA1E Rabbit pAb
Remarks/Alias	BII; CACH6; DEE69; gm139; Cav2.3; EIEE69; CACNL1A6; CACNA1E
Species	Human
GeneID (Human)	777
GeneID	777
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2014-2313 of human CACNA1E (NP_001192222.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15878
Protein Weight	262kDa
Shipping	Ice bag