


Anti-Cacna1H Polyclonal Antibody For If/Icc, Elisa - O95180

Item Number: CM0015466



Introduction

CACNA1H Rabbit pAb is a polyclonal antibody against human CACNA1H, validated for IF/ICC and ELISA, with cross-reactivity to mouse. Designed for research on T-type calcium channel Cav3.2 in ion transport and neurological studies.

[Learn More](#)

Parameter	Value
Product Name	CACNA1H Rabbit pAb
Remarks/Alias	ECA6; EIG6; HALD4; Cav3.2; CACNA1HB; CACNA1H
Species	Human
GeneID (Human)	8912
GeneID	8912
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1890-2030 of human CACNA1H (NP_066921.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	O95180
Protein Weight	259kDa
Shipping	Ice bag