


Anti-Cacna1H Rabbit Polyclonal Antibody For If-P, Elisa - O95180

Item Number: CM0022944



Introduction

Rabbit polyclonal antibody targeting human CACNA1H (Cav3.2), a T-type calcium channel subunit. Validated for IF-P and ELISA, cross-reacts with rat. Ideal for neuroscience, cardiology, and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	CACNA1H Rabbit pAb
Remarks/Alias	ECA6; EIG6; HALD4; Cav3.2; CACNA1HB; CACNA1H
Species	Human
Gene ID (Human)	8912
Gene ID	8912
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 2250-2353 of human CACNA1H (NP_066921.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Rat
SWISS	O95180
Protein Weight	259kDa
Shipping	Ice bag