

Anti-Cacna2D2 Polyclonal Antibody For Wb, Elisa - Q9Ny47

Item Number: CM0017816

Introduction

Rabbit polyclonal antibody targeting human CACNA2D2, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for calcium channel research and related studies.

[Learn More](#)

Parameter	Value
Product Name	CACNA2D2 Rabbit pAb
Remarks/Alias	CASVDD; CACNA2D; CACNA2D2
Species	Human
Gene ID (Human)	9254
Gene ID	9254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-200 of human CACNA2D2 (NP_001167522.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NY47
Protein Weight	130kDa
Shipping	Ice bag