

Anti-Cacna2D3 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Izs8

Item Number: CM0017551

Introduction

Rabbit polyclonal antibody targeting human CACNA2D3 (Voltage-dependent calcium channel subunit alpha-2/delta-3). Validated for WB and ELISA, cross-reacts with human, mouse, rat. For ion channel and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CACNA2D3 Rabbit pAb
Remarks/Alias	HSA272268, CACNA2D3
Species	Human
Gene ID (Human)	55799
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1040 of human CACNA2D3 (NP_060868.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZS8
Protein Weight	123kDa
Shipping	Ice bag