


Anti-Cacnb1 Rabbit Polyclonal Antibody For Wb, Elisa - Q02641

Item Number: CM0025495



Introduction

Rabbit polyclonal antibody targeting CACNB1, a regulatory subunit of L-type calcium channels, validated for WB and ELISA in human, mouse, rat. Ideal for cardiovascular and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CACNB1 Rabbit pAb
Remarks/Alias	CAB1; CCHLB1; CACNLB1; CACNB1
Species	Human
GeneID	782
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 429-598 of human CACNB1 (NP_000714.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02641
Protein Weight	66kDa
Shipping	Ice bag