

Anti-Cacnb3 Polyclonal Antibody For Wb, Elisa - P54284

Item Number: CM0016354



Introduction

Rabbit polyclonal antibody targeting human CACNB3 (CAB3, CACNLB3), a voltage-dependent L-type calcium channel subunit beta-3. Validated for WB and ELISA, reacts with mouse and rat. Ideal for ion channel and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CACNB3 Rabbit pAb
Remarks/Alias	CAB3, CACNLB3, CACNB3
Species	Human
GeneID	784
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 355-484 of human CACNB3 (NP_000716.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P54284
Protein Weight	55kDa
Shipping	Ice bag