

Anti-Cacng1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q06432

Item Number: CM0013526

Introduction

Rabbit polyclonal antibody targeting human CACNG1 (voltage-dependent calcium channel gamma-1 subunit). Validated for WB, IF/ICC, ELISA. Detects human, mouse, rat. Ideal for skeletal muscle L-type calcium channel research.

[Learn More](#)

Parameter	Value
Product Name	CACNG1 Rabbit pAb
Remarks/Alias	CACNLG; CACNG1
Species	Human
GeneID (Human)	786
GeneID	786
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-110 of human CACNG1 (NP_000718.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06432
Protein Weight	25kDa
Shipping	Ice bag