


Anti-Cacybp Polyclonal Antibody For Wb, Elisa - Q9Hb71

Item Number: CM0024398



Introduction

Calcyclin-binding protein (CACYBP) rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Useful for studying ubiquitin-mediated degradation and calcium signaling. Backed by CamelBio's IVD raw-material supply expertise.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CACYBP Rabbit pAb
Remarks/Alias	SIP, GIG5, PNAS-107, S100A6BP, BP
Species	Human
GeneID (Human)	27101
GeneID	27101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of human CACYBP (NP_055227.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HB71
Protein Weight	26kDa
Shipping	Ice bag