


Anti-Cacybp Polyclonal Antibody (Ko Validated) - Q9Hb71

Item Number: CM0027214



Introduction

Rabbit polyclonal antibody against human CACYBP, validated by knockout, for Western blot, IF/ICC, and ELISA. Detects endogenous CACYBP across human, mouse, and rat samples. Recombinant immunogen ensures lot-to-lot consistency.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CACYBP Rabbit pAb
Remarks/Alias	SIP; GIG5; PNAS-107; S100A6BP; BP
Species	Human
GeneID (human)	27101
GeneID	27101
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HB71
Protein Weight	26kDa
Shipping	Ice bag