


Anti-Cad Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P27708

Item Number: CM0029335



Introduction

High-quality rabbit polyclonal antibody targeting human CAD (P27708), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for pyrimidine biosynthesis and metabolic disorder research.

[Learn More](#)

Parameter	Value
Product Name	CAD Rabbit pAb
Remarks/Alias	CDG1Z; DEE50; GATD4; EIEE50; CAD
Species	Human
GeneID (Human)	790
GeneID	790
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1650-1900 of human CAD (NP_004332.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27708
Protein Weight	243kDa
Shipping	Ice bag