


Anti-Cad Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P27708

Item Number: CM0009408



Introduction

Anti-CAD rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Recognizes human CAD, a 243 kDa multifunctional enzyme of the de novo pyrimidine biosynthesis pathway. Ideal for studying pyrimidine metabolism and related disorders.

[Learn More](#)

Parameter	Value
Product Name	CAD Rabbit mAb
Remarks/Alias	CDG1Z; DEE50; GATD4; EIEE50; CAD
Species	Human
GeneID (Human)	790
GeneID	790
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CAD (P27708).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P27708
Protein Weight	243kDa
Shipping	Ice bag