

Anti-Cadherin 16 (Cdh16) Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - O75309

Item Number: CM0027144

Introduction

Rabbit polyclonal antibody against human Cadherin-16 (CDH16, Ksp-cadherin). Validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for kidney-specific cadherin research.

[Learn More](#)

Parameter	Value
Product Name	Cadherin 16 Rabbit pAb
Remarks/Alias	CDH16; Cadherin 16
Species	Human
GenelD	1014
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-230 of human Cadherin 16 (NP_004053.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75309
Protein Weight	90kDa
Shipping	Ice bag