

Anti-Cadherin 16 Monoclonal Antibody For Wb, If-P, Elisa - 075309

Item Number: CM0003602

Introduction

Anti-Cadherin 16 Rabbit Monoclonal Antibody validated for WB, IF-P, ELISA. Targets kidney-specific cadherin, cell adhesion protein. Cross-reacts with Human, Mouse. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	Cadherin 16 Rabbit mAb
Remarks/Alias	Ksp cadherin; Cadherin 16
Species	Human
GenelD (Human)	1014
GenelD	1014
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	075309
Protein Weight	90kDa
Shipping	Ice bag