

Anti-Cadherin-17 (Li Cadherin) Rabbit Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - Q12864

Item Number: CM0011713



Introduction

Rabbit monoclonal antibody targeting human Cadherin-17 (LI-cadherin) validated for WB, IHC-P, IF-P, ELISA, and cross-reactive with mouse and rat. Suitable for cell adhesion and tissue organization studies.

[Learn More](#)

Parameter	Value
Product Name	LI Cadherin/Cadherin-17 Rabbit mAb
Remarks/Alias	HPT1; CDH16; HPT-1; LI Cadherin/Cadherin-17
Species	Human
Gene ID (Human)	1015
Gene ID	1015
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LI Cadherin/Cadherin-17 (Q12864)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12864
Protein Weight	92 kDa
Shipping	Ice bag