

Anti-Cadps Polyclonal Antibody For Wb, Elisa - Q9ULu8

Item Number: CM0013407



Introduction

Rabbit polyclonal antibody targeting human CADPS/CAPS-1, validated for Western blot and ELISA. Cross-reacts with mouse and rat, ideal for neuroscience and endocrine secretion studies involving dense-core vesicle exocytosis.

[Learn More](#)

Parameter	Specification
Product Name	CADPS Rabbit pAb
Remarks/Alias	CAPS; CAPS1; CADPS1; UNC-31; CADPS
Species	Human
Gene ID (Human)	8618
Gene ID	8618
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1020-1100 of human CADPS (NP_899630.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9ULU8
Protein Weight	153kDa
Shipping	Ice bag