

Anti-Calbindin Rabbit Monoclonal Antibody For Wb - P05937

Item Number: CM0004722



Introduction

High-specificity Calbindin Rabbit Monoclonal Antibody for WB, IF-P, ELISA applications. Validated in human, mouse, rat samples. Detects endogenous calbindin at 30 kDa. Ideal for calcium signaling and diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Calbindin Rabbit mAb
Remarks/Alias	CALB; D-28K; Calbindin
Species	Human
GeneID (Human)	793
GeneID	793
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Calbindin (P05937).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05937
Protein Weight	30 kDa
Shipping	Ice bag