


Anti-Calbindin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05937

Item Number: CM0016646



Introduction

Calbindin Rabbit Polyclonal Antibody validated for Western blot (WB), IHC-P, and ELISA. Recognizes human, mouse, and rat Calbindin. Generated with recombinant full-length immunogen. Suitable for calcium signaling, vitamin D, and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Calbindin Rabbit pAb
Remarks/Alias	CALB; D-28K; Calbindin
Species	Human
Gene ID (Human)	793
Gene ID	793
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human Calbindin (NP_004920.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P05937
Protein Weight	30kDa
Shipping	Ice bag