


## Anti-Calbindin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05937

Item Number: CM0015912



### Introduction

Rabbit polyclonal antibody targeting human Calbindin (CALB; D-28K) for Western blot, IHC-paraffin, and ELISA. Cross-reacts with mouse and rat; detects ~30 kDa protein.

[Learn More](#)

| Parameter        | Value                                                                                                                                     |
|------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Calbindin Rabbit pAb                                                                                                                      |
| Remarks/Alias    | CALB; D-28K; Calbindin                                                                                                                    |
| Species          | Human                                                                                                                                     |
| GeneID (Human)   | 793                                                                                                                                       |
| GeneID           | 793                                                                                                                                       |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human Calbindin (NP_004920.1). |
| Source           | Rabbit                                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                                     |
| Application      | WB, IHC-P, ELISA                                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                         |
| SWISS            | P05937                                                                                                                                    |
| Protein Weight   | 30kDa                                                                                                                                     |
| Shipping         | Ice bag                                                                                                                                   |