


Anti-Calca Rabbit Polyclonal Antibody For Ihc-P, Elisa - P06881

Item Number: CM0020097



Introduction

CALCA Rabbit pAb is a rabbit polyclonal antibody targeting human calcitonin gene-related peptide 1 (alpha-CGRP), validated for IHC-P and ELISA, with mouse cross-reactivity. Ideal for vasodilation and neuropeptide research.

[Learn More](#)

Parameter	Value
Product Name	CALCA Rabbit pAb
Remarks/Alias	CT, KC, PCT, CGRP, CALC1, CGRP1, CGRP-I, CGRP-alpha, CALCA
Species	Human
Gene ID (Human)	796
Gene ID	796
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-128 of CALCA (NP_001029125.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P06881
Protein Weight	14kDa
Shipping	Ice bag