

Anti-Calceinurin A Rabbit Monoclonal Antibody For Wb, Elisa - Q08209

Item Number: CM0010764



Introduction

Calcineurin A Rabbit Monoclonal Antibody validated for WB and ELISA; detects human, mouse, rat PPP3CA. Essential for calcium signaling, NFAT-mediated transcription, and T-cell activation research.

[Learn More](#)

Parameter	Value
Product Name	Calcineurin A Rabbit mAb
Remarks/Alias	CALN; CCN1; CNA1; CALNA; DEE91; IECEE; PPP2B; ACCIID; CALNA1; IECEE1; Calcineurin A
Species	Human
Gene ID (Human)	5530
Gene ID	5530
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 422-521 of human Calcineurin A (Q08209).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08209
Protein Weight	59 kDa
Shipping	Ice bag