

Anti-Calcoco1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9P1Z2

Item Number: CM0027366

Introduction

Rabbit polyclonal antibody against CALCOCO1, a transcription coactivator involved in Wnt signaling and nuclear receptor pathways. Validated for WB and ELISA, reacts with human and mouse. Suitable for research on gene regulation.

[Learn More](#)

Parameter	Value
Product Name	CALCOCO1 Rabbit pAb
Remarks/Alias	Cocoa; PP13275; calphoglin; CALCOCO1
Species	Human
GeneID (Human)	57658
GeneID	57658
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P1Z2
Protein Weight	77kDa
Shipping	Ice bag