

Anti-Calcr Polyclonal Antibody For Wb, Elisa - P30988

Item Number: CM0013032



Introduction

CALCR Rabbit pAb is a polyclonal antibody validated for WB and ELISA, targeting human calcitonin receptor (55 kDa). Ideal for research on GPCR signaling, calcium homeostasis, and bone metabolism.

[Learn More](#)

Parameter	Value
Product Name	CALCR Rabbit pAb
Remarks/Alias	CRT; CTR; CT-R; CTR1; CALCR
Species	Human
Gene ID (Human)	799
Gene ID	799
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-148 of human CALCR (NP_001733.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P30988
Protein Weight	55kDa
Shipping	Ice bag