


Anti-Calcrl Rabbit Polyclonal Antibody For Wb, Elisa - Q16602

Item Number: CM0029858



Introduction

High-quality anti-CALCRL rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, rat CALCRL (CRLR; CGRPR), a key receptor for calcitonin gene-related peptide and adrenomedullin. Ideal for cardiovascular and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	CALCRL Rabbit pAb
Remarks/Alias	CRLR; CGRPR; LMPHM8; CALCRL
Species	Human
Gene ID (Human)	10203
Gene ID	10203
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 346-445 of human CALCRL (NP_005786.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16602
Protein Weight	53kDa
Shipping	Ice bag