

Anti-Caldesmon/Cdm Rabbit Monoclonal Antibody For Ihc-P, Elisa - Q05682

Item Number: CM0012650

Introduction

High-quality rabbit monoclonal antibody against Caldesmon/CDM, validated for IHC-P and ELISA. Detects human Caldesmon, an actin-binding protein regulating smooth muscle and nonmuscle cell contractility, plus mitosis and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Caldesmon/CDM Rabbit mAb
Remarks/Alias	CDM; HCAD; LCAD; h-CD; H-CAD; L-CAD; NAG22
Species	Human
GeneID (Human)	800
GeneID	800
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 208-436 of human Caldesmon/CDM (NP_149129.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q05682
Protein Weight	120-150kDa
Shipping	Ice bag