


Anti-Calhm1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8lu99

Item Number: CM0025147



Introduction

High-quality rabbit polyclonal antibody targeting human CALHM1 (FAM26C), validated for WB, IHC-P, and ELISA. Cross-reacts with human and rat. Suitable for studying taste perception, calcium homeostasis, and neuronal signaling.

[Learn More](#)

Parameter	Value
Product Name	CALHM1 Rabbit pAb
Remarks/Alias	FAM26C; CALHM1
Species	Human
GeneID (Human)	255022
GeneID	255022
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8IU99
Protein Weight	38kDa
Shipping	Ice bag