


Anti-Calhm2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ha72

Item Number: CM0014652



Introduction

Anti-CALHM2 Rabbit Polyclonal Antibody validated for Western blot and ELISA. Detects human and rat CALHM2, a calcium homeostasis modulator channel involved in ATP release, neuroinflammation, and synaptic transmission.

[Learn More](#)

Parameter	Value
Product Name	CALHM2 Rabbit pAb
Remarks/Alias	FAM26B; CALHM2
Species	Human
Gene ID (Human)	51063
Gene ID	51063
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human CALHM2 (NP_057000.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9HA72
Protein Weight	36kDa
Shipping	Ice bag