


Anti-Calm2 Polyclonal Antibody For Wb And Elisa - P0Dp24

Item Number: CM0028309



Introduction

Rabbit polyclonal antibody targeting human CALM2 calmodulin-2, validated for WB and ELISA with cross-reactivity to mouse. Essential for calcium signaling studies, centrosome biology, and cardiac research.

[Learn More](#)

Parameter	Value
Product Name	CALM2 Rabbit pAb
Remarks/Alias	caM; CALM; CAM1; CAM3; CAMC; PHKD; CAMII; LQT15; PHKD2; CALML2; CAMIII; CALM2
Species	Human
GeneID (human)	805
GeneID	805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-149 of human CALM2 (NP_001734.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P0DP24
Protein Weight	17kDa
Shipping	Ice bag