

Anti-Calm3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P62158

Item Number: CM0012903

Introduction

Rabbit polyclonal antibody against human CALM3 (calmodulin 3), validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. 16kDa protein. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	CALM3 Rabbit pAb
Remarks/Alias	CaM; CALM; CAM1; CAM2; CAMB; PHKD; CPVT6; LQT16; PHKD3; CaMIII; HEL-S-72; CALM3
Species	Human
GeneID (human)	808
GeneID	808
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-140 of human CALM3 (NP_005175.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62158
Protein Weight	16kDa
Shipping	Ice bag