

Anti-Calmegin (Clgn) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - O14967

Item Number: CM0002108

Introduction

High-specificity rabbit monoclonal antibody against human Calmegin (CLGN), suitable for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for reproductive biology studies on spermatogenesis and male fertility.

[Learn More](#)

Parameter	Value
Product Name	Calmegin Rabbit mAb
Remarks/Alias	CLGN; Calmegin
Species	Human
GeneID (Human)	1047
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-350 of human Calmegin (O14967)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14967
Protein Weight	70kDa
Shipping	Ice bag