

Anti-Calmodulin 1/2/3 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P0Dp23 / P0Dp24 / P0Dp25

Item Number: CM0013335



Introduction

Rabbit polyclonal antibody targeting human Calmodulin 1/2/3, suitable for IHC-P, IF/ICC, and ELISA. Detects all three calmodulin isoforms in human, mouse, and rat. Ideal for calcium signaling and channel regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Calmodulin 1/2/3 Rabbit pAb
Remarks/Alias	CALM1; CALML2; CAMI; CPVT4; DD132; LQT14; PHKD; caM; Calmodulin 1/2/3
Species	Human
Gene ID (Human)	801/805/808
Gene ID	801/805/808
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 70-149 of human Calmodulin 1/2/3 (NP_008819.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DP23/P0DP24/P0DP25
Protein Weight	16 kDa
Shipping	Ice bag