


Anti-Calmodulin 1/2/3 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P0Dp23 / P0Dp24 / P0Dp25

Item Number: CM0012234



Introduction

Rabbit monoclonal antibody targeting Calmodulin 1/2/3, validated for WB, IF-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for calcium signaling and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Calmodulin 1/2/3 Rabbit mAb
Remarks/Alias	CALML2; CAMI; CPVT4; DD132; LQT14; PHKD; caM; Calmodulin 1/2/3
Species	Human
Gene ID (Human)	801/805/808
Gene ID	801/805/808
Immunogen	Synthetic Peptide. A synthetic peptide corresponding to a sequence within amino acids 70-149 of human Calmodulin (NP_008819.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DP23/P0DP24/P0DP25
Protein Weight	16 kDa
Shipping	Ice bag