

Anti-Calmodulin 1/2/3 Rabbit Polyclonal Antibody For Wb, Elisa - P0Dp23/P0Dp24/P0Dp25

Item Number: CM0016823



Introduction

Calmodulin 1/2/3 rabbit polyclonal antibody validated for Western blot (WB) and ELISA applications. Targets human, mouse, rat calmodulin, a key calcium-binding signal transducer. Recombinant antibody from rabbit for research use.

[Learn More](#)

Parameter	Value
Product Name	Calmodulin 1/2/3 Rabbit pAb
Remarks/Alias	CALM1; CALML2; CAMI; CPVT4; DD132; LQT14; PHKD; caM; Calmodulin 1/2/3
Species	Human
GeneID (Human)	801/805/808
GeneID	801/805/808
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 70-149 of human Calmodulin 1/2/3 (NP_008819.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DP23/P0DP24/P0DP25
Protein Weight	16kDa
Shipping	Ice bag