

Anti-Calmodulin (Calm3) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62158

Item Number: CM0015505

Introduction

Rabbit polyclonal antibody against human calmodulin (CALM3) validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. Ideal for calcium signaling and cell regulation studies.

[Learn More](#)

Parameter	Value
Product Name	CALM3 Rabbit pAb
Remarks/Alias	CaM; CALM; CAM1; CAM2; CAMB; PHKD; CPVT6; LQT16; PHKD3; CaMIII; HEL-S-72; CALM3
Species	Human
GeneID (human)	808
GeneID	808
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-149 of human CALM3 (NP_001734.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P62158
Protein Weight	17kDa
Shipping	Ice bag