

Anti-Calnexin Monoclonal Antibody For Wb, Elisa - P27824

Item Number: CM0010513

Introduction

Calnexin Rabbit mAb for Western blot and ELISA. Reacts with human samples. Targets 68kDa calnexin protein, a key ER chaperone involved in glycoprotein quality control and protein folding. For cellular and molecular biology research.

[Learn More](#)

Parameter	Value
Product Name	Calnexin Rabbit mAb
Remarks/Alias	CNX; P90; IP90; Calnexin
Species	Human
Gene ID (Human)	821
Gene ID	821
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 493-592 of human Calnexin (P27824).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P27824
Protein Weight	68kDa
Shipping	Ice bag