

## Anti-Calnexin Polyclonal Antibody For Wb, If, Ihc, Elisa - P27824

Item Number: CM0018069

### Introduction

Rabbit polyclonal antibody targeting human Calnexin (CNX, IP90, P90), an ER chaperone. Validated for WB, IF/ICC, IHC-P, and ELISA. Detects Calnexin in human, mouse, and rat samples. Ideal for protein folding and quality control research.

[Learn More](#)

| Parameter        | Value                                                                                                      |
|------------------|------------------------------------------------------------------------------------------------------------|
| Product Name     | Calnexin Rabbit pAb                                                                                        |
| Remarks/Alias    | CNX; P90; IP90; Calnexin                                                                                   |
| Species          | Human                                                                                                      |
| GeneID (Human)   | 821                                                                                                        |
| GeneID           | 821                                                                                                        |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 500-592 of human Calnexin (NP_001737.1) |
| Source           | Rabbit                                                                                                     |
| Category         | Polyclonal Antibodies                                                                                      |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                          |
| SWISS            | P27824                                                                                                     |
| Protein Weight   | 68 kDa                                                                                                     |
| Shipping         | Ice bag                                                                                                    |