


Anti-Calpain 1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07384

Item Number: CM0018254



Introduction

Rabbit polyclonal antibody against human calpain-1 (CAPN1), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studies on cytoskeletal remodeling and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	Calpain 1 Rabbit pAb
Remarks/Alias	CANP; muCL; CANP1; SPG76; CANPL1; muCANP; Calpain 1
Species	Human
GeneID (Human)	823
GeneID	823
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 355-530 of human Calpain 1 (NP_005177.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07384
Protein Weight	82kDa
Shipping	Ice bag