

Anti-Calpain 1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P07384

Item Number: CM0012767



Introduction

Rabbit monoclonal antibody against human Calpain 1 (muCANP), validated for WB, IF/ICC, ELISA in human, mouse, rat. Suitable for apoptosis and cytoskeletal remodeling studies.

[Learn More](#)

Parameter	Value
Product Name	Calpain 1 Rabbit mAb
Remarks/Alias	CANP; muCL; CANP1; SPG76; CANPL1; muCANP; Calpain 1
Species	Human
GeneID (Human)	823
GeneID	823
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human Calpain 1 (P07384).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07384
Protein Weight	82 kDa
Shipping	Ice bag