


Anti-Calpain 1 Rabbit Polyclonal Antibody For Wb, Elisa - P07384

Item Number: CM0021313



Introduction

High-quality rabbit polyclonal antibody targeting human Calpain 1 (muCANP/CANP1). Validated for WB and ELISA, with a protein weight of 82kDa. Ideal for investigating calcium-regulated proteolysis and neurodegenerative pathways.

[Learn More](#)

Parameter	Value
Product Name	Calpain 1 Rabbit pAb
Remarks/Alias	CANP; muCL; CANP1; SPG76; CANPL1; muCANP; Calpain 1
Species	Human
Gene ID (Human)	823
Gene ID	823
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-330 of human Calpain 1 (NP_005177.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P07384
Protein Weight	82kDa
Shipping	Ice bag