

Anti-Calpain 2 Rabbit Monoclonal Antibody For Wb, Elisa - P17655

Item Number: CM0010151



Introduction

Rabbit monoclonal antibody against human Calpain-2 catalytic subunit (CANP2, mCANP). Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for cytoskeletal remodeling and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Calpain 2 Rabbit mAb
Remarks/Alias	CANP2, mCANP, CANPL2, CANPml, Calpain 2
Species	Human
Gene ID (Human)	824
Gene ID	824
Immunogen	A synthetic peptide corresponding to amino acids 600-700 of human Calpain-2 (P17655)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17655
Protein Weight	80kDa
Shipping	Ice bag