

Anti-Calpain 2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17655

Item Number: CM0027539



Introduction

Rabbit polyclonal antibody against human Calpain 2. Suitable for WB, IHC-P, and ELISA; cross-reacts with human, mouse, and rat. Ideal for studying calpain-mediated proteolysis and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	Calpain 2 Rabbit pAb
Remarks/Alias	CANP2; mCANP; CANPL2; CANPml; Calpain 2
Species	Human
GeneID (Human)	824
GeneID	824
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 520-700 of human Calpain 2 (NP_001739.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB,IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17655
Protein Weight	80kDa
Shipping	Ice bag