


Anti-Calpain Small Subunit 1 (Capns1) Rabbit Monoclonal Antibody For Wb, Elisa - P04632

Item Number: CM0011221



Introduction

Rabbit monoclonal antibody against human Calpain small subunit 1 (CAPNS1, CAPN4). Suitable for Western blot and ELISA; cross-reacts with mouse and rat. For research on calcium-regulated proteolysis and cytoskeletal dynamics.

[Learn More](#)

Parameter	Value
Product Name	Calpain small subunit 1 Rabbit mAb
Remarks/Alias	CANP; CDPS; CSS1; CANPS; CAPN4; CALPAIN4; Calpain small subunit 1
Species	Human
GeneID (Human)	826
GeneID	826
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Calpain small subunit 1 (P04632)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04632
Protein Weight	28kDa
Shipping	Ice bag