


Anti-Calponin 1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P51911

Item Number: CM0005579



Introduction

Rabbit monoclonal antibody against human Calponin 1 (CNN1, SMCC). Validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for smooth muscle contraction and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	Calponin 1 Rabbit mAb
Remarks/Alias	SMCC; Sm-Calp; HEL-S-14
Species	Human
Gene ID (Human)	1264
Immunogen	Synthetic peptide (aa 198-297 of human Calponin 1, NP_001290.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51911
Protein Weight	33kDa
Shipping	Ice bag