


Anti-Calponin-1 Rabbit Monoclonal Antibody For Wb, Ihc, If, Elisa - P51911

Item Number: CM0010147



Introduction

Recombinant rabbit monoclonal antibody targeting human Calponin-1 (CNN1, SMCC). Validated in WB, IHC-P, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for smooth muscle contraction studies.

[Learn More](#)

Parameter	Value
Product Name	Calponin Rabbit mAb
Remarks/Alias	SMCC; Sm-Calp; HEL-S-14; Calponin
Species	Human
Gene ID (Human)	1264
Gene ID	1264
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 198-297 of human Calponin (P51911).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51911
Protein Weight	33kDa
Shipping	Ice bag